


EASYTRACK® WALL MOUNT ADAPTOR

MODEL #700.10860

INSTALLATION AND PRODUCT CARE INSTRUCTIONS



IMPORTANT SAFETY DIRECTIONS

WARNING: Failure to observe Safety Directions listed below could result in serious injury.

Before you install the Easytrack system:

- ◆ Carefully inspect product upon arrival. If any component is missing (see Figure 8) or appears to be damaged, do not install the Easytrack system. Arrange for repair/replacement of any faulty parts by contacting your local distributor for further instructions.
- ◆ If you are unsure of the strength or stability of your wall, do not install the Wall Mount Adaptor. Consult your local distributor.
- ◆ Do not install on a wall with water damage.
- ◆ Do not drop product, because it may cause breakage.
- ◆ Do not completely submerge the Easytrack system in water.
- ◆ Do not attempt to over-extend the rail past the locking point.

WARNING:

- ◆ Do not use the Easytrack system as a swing.
- ◆ Refrain from swinging patient sideways of rail.

GETTING STARTED

Upon receipt of your Easytrack Wall Mount Adaptor:

WARNING: Read the "Requirements" and "Important Safety Directions" sections before attempting to assemble the Wall Mount Adaptor.

- ◆ Survey the room where you plan to install the Easytrack Wall Mount Adaptor.
- ◆ Make sure the room meets all the specifications outlined in the section "Requirements".
- ◆ Use the measurements below to determine the best place to put the Easytrack system :
 - ◆ Over a bed – approx. 915 – 1020 mm/36 - 40 in from the head of the bed (see Figure 1);
 - ◆ Over a tub – approx. 400 – 460 mm/16 – 18 in from the end of the tub (see Figure 1);
 - ◆ Allow 915 mm/ 36 in of open area for transfer working space (see Figure 2);
 - ◆ Allow a minimum of 1340 mm/53 in from the bottom of the rail to the top of any object to be transferred over (see Figure 3).




Figure 1

Figure 2




Figure 3

- ◆ Transfers should only take place under the Easytrack rail (see **Figure 4** and **Figure 5**).




Figure 4




Figure 5

- ◆ The Wall Mount Adaptor may be installed on either a concrete wall or wood-framed wall. For wood, the Wall Mount Adaptor must be attached to a minimum of three (3) 2 X 3 studs. DO NOT attempt to install the wall mount on less than three (3) studs (see **Figure 6**).




Figure 6

- The Wall Mount Adaptor must be positioned on the wall to ensure that the angle between the wall and the rail will measure more than 30° (see Figure 7).




Figure 7

WHAT’S INCLUDED (WALL MOUNT ADAPTOR)

Item	Description	Part number	Qty
1	Wall mount plate	211.10850	1
2	Pivoting rail attachment	700.10870	1
Includes:			
3	Round blue cap	200.10720	1
4	Pivot wall mount	N/A	1
5	M10 x 25 screw	000.04090	4
6	Lockwasher M10	000.04425	4
7	Fixation wall mount	N/A	1
8	Back plate	N/A	1




Figure 8

ASSEMBLY

ASSEMBLY FOR WALL BRACKET EASYTRACK SYSTEM

WARNING: Read the section "Requirements" and the warnings before attempting to assemble the Wall Mount Adaptor system.

REQUIRED ELEMENTS

- Step ladder
- Tape measure
- Six (6) concrete 3/8" anchors or six (6) 3/8" lag bolts
- Drill
- Wrench with a 9/16" socket
- Stud finder
- Someone to assist you

- 1) Refer to Figures 9 and 10 below for the distance required between the ceiling and where the top of the Wall Mount Adaptor must be installed. Please note that when the wall bracket has the pivoting rail of a 3-post system attached, the mounting point must be lower than a standard rail.




Figure 9

Figure 10

- 2) Install the wall plate on the wall. The wall plate must be installed on a concrete wall with six (6) 3/8 in concrete anchors or on three (3) 2 x 3 (minimum size) wooded studs with two (2) 3/8" dia. (M10) and 4 in long lag bolts on each stud (see Figure 11).




Figure 11

NOTE: Make sure the rail is levelled before securing it to the wall.

Wood studs: Use a stud finder to find the center of each stud. It is necessary to pre-drill holes into each stud before screwing the lag bolt to prevent the wood from cracking.

Follow drill size recommendation for lag bolt:

LAG BOLT	M10 (3/8")
PRE-DRILL SIZE	6.5 mm (1/4")

WARNING: Pre-drilled holes must be in the center of the stud.
If holes are not properly centered, serious damage or injury may occur.

WARNING: DO NOT install steel studs, since reinforcement may be required.
Call your ArjoHuntleigh dealer for more information.

WARNING: Avoid drilling close to any electrical wiring or plumbing. Please make sure you know what is behind the wall you are drilling in before doing any modifications to it. Failure to do so may cause serious injury.

- 3) Insert the pivoting rail attachment into the wall plate. Make sure that the bracket is placed properly (see Figure 12).




Figure 12

- 4) See the section "Getting Started" for suggested positions. Slide the pivoting rail attachment into its final position. For the 3-post assembly, the bracket must be placed in a way that the rail will not, when turning, generate an angle beyond the maximum allowed (see Figure 13).




Figure 13

- 5) Tighten the four (4) bolts on pivoting rail attachment. The rail can now be attached to the Wall Mount Adaptor.

SPECIFICATIONS

EASYTRACK WALL MOUNT ADAPTOR SPECIFICATIONS

SPECIFICATIONS	
Product code	700.10860
Capacity	200 kg (440 lb)
Material	Anodised aluminium and steel
Weight	12 kg (27 lb)

WALL MOUNT ADAPTOR DIMENSIONS




Figure 14

MAINTENANCE

WARNING: Failure to follow maintenance instructions below when indicated could result in injury or damage.

WARNING: Periodic inspection should be performed by a person who is suitable and properly qualified and well acquainted with the design, use and care of the Easytrack system. Contact your local dealer/distributor.

WARNING: Inspection should be carried out routinely (at least as often as specified), and as necessary. All defects and damages that have lead to corrective actions should be noted, signed and dated by the inspector. Defects and resulting actions should be reported to your supplier.

Before each time you use the Easytrack Wall Mount Adaptor:

1. Check that the rail is secured in place. Press the rail upward to make sure that the rail is locked onto the post.
2. Check if there is any appearance of damage. IF THERE IS ANY APPEARANCE OF DAMAGE, DO NOT USE THE EASYTRACK SYSTEM.

Every month the following should be done :

1. Clean inside of rail with a damp cloth to ensure there is no dust or grease build-up inside of rail.
2. Check to ensure you still have all your original parts. If any part is missing, it should be replaced BEFORE USING THE EASYTRACK WALL MOUNT ADAPTOR AGAIN.
3. Check for the appearance of rust. If rust is starting to appear, simply scrub the piece clean. Examine all joints and pieces for wear and fatigue. Replace any part that looks damaged or worn.
4. Keep a log of any repair/replacement/inspection.