



Etac Swift Mobile Tilt

With a unique tilt in space function.

Etac Swift Mobile Tilt is available in two models. Both have manually operated tilt functions but with different weight capacity. They all have the same characteristic design and functions as Etac Swift Mobile, page 28.

Tilting operated from the side

Etac Swift Mobile Tilt has the unique function which enables tilt manoeuvring from the side. This gives the best conditions for an ergonomically correct work position for the carer. It also enables eye contact and communication with the care recipient during the tilting process.

Assisted movement back to an upright position

The gas pistons work to facilitate bringing the seat into an upright position. When tilting, the user's body weight will act as counterweight and make the movement smooth. This reduces the risk of carers suffering work related injuries. The seat height can be set in various fixed positions and it's done without using tools.

A tilted position increases stability and facilitates in care situations.



Description	Item no.	Materials
Swift Mobile Tilt	80209410	Gas pistons: lift capacity 130N
Swift Mobile Tilt with bucket holder	80209411	Other material: see "Etac Swift Mobile shower commode" on page 28
Swift Mobile Tilt with pan holder	80209412	

Standard: head support, back support, arm supports, seat, foot supports, heels traps, and 4 lockable castors.






Dimensions

Adjustable tilt in space: 0° to -35° can be set to +5° to -30°


Weight

19.7 kg (43.4 lbs)
With power unit: 23.6 kg (52 lbs)

Maintenance



Clean the product with a solvent-free cleaning agent (pH level 5–9), or with a 70% disinfectant solution. Can be decontaminated at max 85 °C / 185 °F if so required.




Features



Unique tilt bar

The bar for the tilt mechanism is easy to reach from either side.



Head Support

The adjustability ensures a comfortable and restful support.



Close contact during the tilting process

The unique tilt function enables an ergonomically correct working position for the carer who is able to maintain eye contact with the care recipient.

