



**Wave** (Z220 with  
tub transfer base Z228)



# Wave features



**Generously sized**  
The wide seat comfortably accommodates large clients. The seat width on the large size is 17½".



**One-hand adjustment**  
Adjust backrest angle with one hand, keeping the other hand on the client for safety and reassurance.



**Full leg support**  
The calf rest provides full support for the lower legs. Adjust it independently of the seat by pushing the gray buttons and rotating up or down. It must be installed at the time of purchase.



**Typical bathtub height range**  
—— 18"  
—— 14"





**Tub stand**



**Shower stand**



**Tub transfer base**



**Highly adjustable**

Backrest angle, seat angle and calf rest angle adjust individually.



**More seat heights**

For an intermediate seat height, put the rear legs up and front legs down. To get the chair low in the tub, put the front and back legs up. Fold the chair completely flat for storage.



# Components

## Bath chair size and color



## Choice of chest strap



### Without lateral positioning

The simple chest strap safely positions the user in the Wave bath chair. It attaches with hook and loop where needed along the full height of the back and buckles securely behind.



### With lateral positioning

We designed this chest strap with adjustable lateral straps to give additional lateral support, helping users maintain better trunk alignment in the bath chair. The soft fabric laterals will not get in the way during transfers.







### Head blocks

Head blocks provide lateral support and positioning for the user's head.



### Leg straps

Use the fully adjustable leg straps to position your client's hips in abduction (or just to provide additional security).



### Calf rest

The calf rest supports the lower legs. It is strongly recommended when the bath chair is used on one of the stands.

**Note:** The calf rest must be chosen at the time of the original purchase. It cannot be added later.



### Tub stand

The tub stand raises the bath chair either 5½" or 11½" higher in the tub depending on the position of the tub stand legs.



### Shower stand

The shower stand lets you roll your client into the shower stall at a convenient height for cleaning. Locking 4" casters provide smooth, easy maneuverability.



## Components (continued)

### Tub transfer base

Bathing equipment for people with disabilities tends to be big and bulky, taking up far too much space in bathrooms shared by the whole family. Using Rifton's ingenious new tub transfer base, you can perform healthy transfers without crowding your bathroom with permanently mounted equipment. When bath time is over, the base can be folded up and hung compactly on the wall.



### Unique two-stage movement



The tub transfer base starts in the position farthest out of the tub. It is then rolled back over the tub, and the brakes are set.



As the upper portion with the bath chair is rotated 90° it moves approximately 12" back over the tub in one motion.



The legs of the tub transfer base can be adjusted to fit bathtubs up to 24" high.





**1. With the client secured in the chair,** push the transfer base to the back of the tub as far as it will go and lock the wheels.



**2. Press either of the gray pivot latches** and rotate the bath chair and client over the tub.



**3. As the bath chair rotates,** it continues to move back over the tub. Gently lift the client's legs to clear the edge of the tub.



**4. In the final position** the bath chair is over the tub and can be adjusted as desired for showering.

## Free up the tub.

**Rifton's folding,** storable tub transfer base makes it possible for everyone in the family to still use the tub.

When you are finished, remove the base from the tub, fold the legs and hang it on the wall hook provided.

