



A



A

Bariatric Side Walker

Platform Walker/Crutch Attachment

Item #	Product Description	UOM
10242	Adult (Hemi-Walker)	1/ea

- Designed for individuals with the use of only one hand or arm.
- Lighter than a walker and more stable than a cane.
- Folds easily with one hand. (Figure A)
- Ideal for users who do not require a walker but need a wider base for support.
- 500 lb. Weight Limit.

specifications	
Overall Dimensions	18.5" (W) x 16" (D)
Height	29.5" - 36"
Weight	4 lbs.
Carton Shipping Weight	5 lbs.
Weight Limit	500 lbs.
Warranty	Limited Lifetime

Item #	Product Description	For Use With	Product Height	UOM
10105-1		Adult and Junior Walkers and Aluminum Crutches	37" - 53"	1/ea
10105-2		Adult and Junior Walkers and Aluminum Crutches	37" - 53"	2/cs
10116	Universal Platform Walker/Crutch Attachment	All Drive walkers and Invacare® Paddle Release Walkers	36" - 55"	1/pr

- Designed for those who are unable to grip the walker.
- Cradles forearm with soft vinyl padding for comfort.
- Large top bracket allows the platform attachment to be placed over the hand grip of the walker.
- Adjustable strap securely holds arm in place.
- 10105CB (Crutch Attachment) sold separately.

specifications	
Overall Width (10105 10116)	5.75" 6"
Overall Depth	18"
Weight (10105 10116)	2 7 lbs.
Carton Shipping Weight (1010501 10105-2 10116)	3 6 10 lbs.
Weight Limit	300 lbs.
Warranty	Limited Lifetime