



## Airgo Comfort-Plus Adjustable Quad Canes, Large Base

HCPCS:




FIGURE A

- For people who need extra help with balance, our quad cane offers added support in a lightweight design
- Built for longevity from durable anodized aluminum
- The K-shaped rotating base swivels and locks into position to accommodate right-handed or left-handed users and is available in small or large size for greater stability
- Ergonomic cushion-top handle conforms to your hand and helps reduce impact shock
- Designed to fit more comfortably in your hand and is ideal for people with arthritis
- Cane strap with reflective patch for added night time safety
- 11 individual height settings offer an unparalleled range of adjustments that fits most people 5' - 6'5"
- Free Standing – Conveniently stands on its own
- AADL# C086

PRODUCTS	Description
770-842	Airgo Comfort-Plus Adjustable Quad Canes, Large Base, EA
770-844	Airgo Comfort-Plus Adjustable Quad Canes, Large Base, EA
770-846	Airgo Comfort-Plus Adjustable Quad Canes, Large Base, EA

SPECIFICATIONS	Description
Base Size	7 in x 11 in (17.8 x 27.9 cm)
Seat Size	20" x 12" (51 x 31 cm)
Size	31.5 in x 15.75 in (80 x 40 cm)
Warranty	Limited Lifetime
Weight	4.4 lbs. (2 kg)
Weight Capacity	300 lbs
Weight Capacity (13046   13046-HD)	300 lbs (136 Kg)