



## Airgo Comfort-Plus Aluminum Cane, Derby Handle

HCPCS:



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E

- For people who could use a hand with balance while walking, our Airgo derby handle cane offers the perfect amount of stability and support
- The Comfort-Plus handle with a cushioned grip helps absorb shocks
- It is ideal for those with arthritis
- NEW cane handle grips to the wall and won't fall
- The 12 push-button height settings ensure this cane can be adjusted to fit you properly
- Built for longevity from ultra-durable lightweight anodized aluminum
- Includes a cane strap with reflective patch for added night time safety
- The treaded rubber Flex-Grip tip provides superior traction
- Fits most people 5' ä 6'5"

### PRODUCTS

#### Description

730-462	Airgo Comfort-Plus Aluminum Cane, Derby Handle, EA
730-464	Airgo Comfort-Plus Aluminum Cane, Derby Handle, EA
730-466	Airgo Comfort-Plus Aluminum Cane, Derby Handle, EA
730-468	Airgo Comfort-Plus Aluminum Cane, Derby Handle, EA
730-470	Airgo Comfort-Plus Aluminum Cane, Derby Handle, EA
730-472	Airgo Comfort-Plus Aluminum Cane, Derby Handle, EA

### SPECIFICATIONS

#### Description

Warranty	Limited Lifetime
Weight Capacity	300 lbs
Weight Capacity (13046   13046-HD)	300 lbs (136 Kg)