

•

•

٠

•

Gel Grip, Aluminum Offset Canes, Height Adjustable

Item #	Color	UOM
RTL10372BK	Black	1/ea
RTL10372RC	Red Crackle	1/ea
RTL10372BC	Blue Crackle	1/ea

- Ergonomically designed handle with soft Gel Grip provides comfort and security.
- Comes standard with premium Gel Grip and wrist strap.
- Manufactured with sturdy, extruded aluminum tubing.
- Easy-to-use, one-button height adjustment with locking ring. Locking ring prevents rattling.
- Available in Blue Crackle (main image), Red Crackle (Figure A) or Black (Figure B)
- Handle height adjusts from 30" to 39".
- 300 lb. Weight Capacity.

Aluminum Ortho K-Grip Canes

ltem #	Color	UOM
10350-6	Bronze	6/cs
10350-1	Bronze	1/cs

- Angle-Adjustable Ergonomically designed contoured grip adjusts to fit the user's wrist and hand. This can reduce stress on the user's hands and joints. (Figure A - C)
- Hand Grip can be rotated for left or right-hand use.
- Manufactured with sturdy, extruded aluminum tubing.
- Easy-to-use, one-button height adjustment with
- locking ring. Locking ring prevents rattling. (Figure D) Handle height adjusts from 30" to 39".
- Comes standard with convenient wrist strap.
- 300 lb. Weight Capacity.

specifications		specifications	
Handle Height	30" - 39"	Handle Height	
Weight Limit	300 lbs.	Weight Limit	
Warranty	Limited Lifetime	Warranty	

30" - 39" 300 lbs. Limited Lifetime