



Quad Canes, Large and Small Base, with Silver Vein Finish

Item #	Product Description	UOM
10385-4	Large Base	4/cs
10384-1	Small Base	1/cs
10384-4	Small Base	4/cs

- Tab-Loc silencer provides added safety to ensure a locked height adjustment. (Figure A)
- Features a padded grip. The flat side of the grip faces up to reduce stress that the user's hand may experience.
- Four point base, combined with offset handle provides additional stability and support.
- Attractive silver vein finish on base. (Figure B)
- Ergonomically designed padded handle.
- Manufactured with an aluminum shaft and a steel base to provide the perfect balance of safety and durability.
- Handle height adjusts from 29" to 38".
- 300 lb. Weight Capacity.

specifications

Handle Height	29" - 38"
Weight Limit	300 lbs.
Warranty	Limited Lifetime

Quad Canes, Large Base

Item #	Product Description	Color	UOM
10300-4	Vinyl Contoured Grip	Chrome	4/cs
10300F-4	Foam Rubber Grip	Chrome	4/cs
10308-4	Vinyl Contoured Grip	Black	4/cs
10309-4	Vinyl Contoured Grip	Bronze	4/cs

- Four point base, combined with offset handle provides additional stability and support.
- Easy-to-use, one-button height adjustment with locking bolt, allows quick, safe fitting.
- Available in Black (main image), Chrome (Figure A), Bronze (Figure B) or Floral. (Figure C)
- No assembly required.
- Handle height adjusts from 30" to 39".
- Handle available in foam rubber (main image) or vinyl. (Figure D)
- 300 lb. Weight Capacity.

specifications

Handle Height	30" - 39"
Weight Limit	300 lbs.
Warranty	Limited Lifetime