



The Invacare Deluxe Three-Position Recliner was designed for maximum comfort, from special padding technology and lumbar support to an articulating headrest and convenient tray which easily folds out of the way. It offers full recline, leg-ottoman elevated and upright positions.

Footrests fold under chair, allowing users to rise to a standing position more easily

Blow-molded side panels are easy to remove and clean

Tabs on convenient overlap tray allow the tray to move freely for ease of use

Articulating headrest makes it easier for user to see visitors

Flame-retardant upholstery is also resistant to fluids for peace-of-mind (CA117)

Color:	Rosewood, Blueberry, Green
Dimensions:	47" H x 28" W x 41" D Seat Dimensions: 19" W x 20" D Tray Dimensions: 25.5" L x 15" W
Seat To Floor Height:	21"
Seat To Back Height:	34"
Seat Angle:	Seat Angle Range: 11-35 degrees
Back Angle:	Back Angle Range: 38-61 degrees
Product Weight Capacity:	250 lb.
Warranty:	Limited Warranty - 3 Year - Chair, 1 Year - Upholstery