




Etac Beauty body care

Using the Beauty products not only makes washing the body easier, they also offer a smooth massage for the skin and hair scalp.

Etac Beauty body and back washer

Beauty body and back washers are designed to facilitate reaching and washing. The wash cloths do not absorb water which make the products light and balanced even when wet. The handles are oval and have a surface which enables a secure grip. The wash cloths are machine washable and replaceable.



Etac Beauty Body washer has a user friendly curved design.



Etac Beauty Back washer is designed to facilitate a two hand grip.



Etac Beauty hair washer

Beauty Hair washer substitutes the fingers to apply shampoo and massage the scalp. The shaft and handle are of the same design as the body washer, making it easier to reach the head without strain on hands, shoulders or arms.



Description

| | |
|-------------------------------------|----------|
| Beauty body washer | 80210076 |
| body washer spare wash cloth 2/pack | 80210103 |
| Beauty back washer | 80210080 |
| back washer spare wash cloth | 80210104 |
| Beauty hair washer | 80210078 |

Item no.

Weight

| |
|-------------------------------------|
| Beauty body washer: 112 g (4 oz.) |
| Beauty back washer: 203 g (7.2 oz.) |
| Beauty hair washer: 112 g (4 oz.) |

Materials

Polypropylene and TPE
Wash cloth: Cotton with polythene inner lining

Design

Ergonomidesign
The Beauty product designs are registered



Etac Beauty comb and hairbrush

Beauty Comb and Hairbrush come in two lengths. The shorter version is recommended when the user can reach the ear. The longer version is suitable for users who can only reach the shoulder. Both the comb and brush are curved to ensure the best comfort and control when combing/brushing the hair.



Etac Beauty multipurpose grip

Beauty multipurpose grip rests comfortably in the hand and offers a variety of uses. Put a nail file in the multipurpose grip for an easy pedicure. The wider end of the grip may be used for a tooth brush or shaver.

The Etac Beauty multipurpose grip has a soft smooth surface.



Description

| | |
|--------------------------|----------|
| Beauty hairbrush | 80210072 |
| Beauty hairbrush, long | 80210074 |
| Beauty comb | 80210068 |
| Beauty comb, long | 80210070 |
| Beauty multipurpose grip | 80210082 |

Length

Beauty hairbrush: 29/37 cm (11½/ 14½")
 Beauty comb: 30/37 cm (11¾/ 14½")
 Beauty multipurpose grip: 11 cm (4¼")

Item no.

Weight

Beauty hairbrush 102/112 g (3.6/ 4 oz.)
 Beauty comb 80/90 g (2.8/ 3.2 lbs)
 Beauty multipurpose grip: 45 g (1.6 oz.)

Materials

Polypropylene and TPE
 Nail file: Emery board

Design

Ergonomidesign
 The Beauty product designs are registered

