



## Pivot & Tilt Overbed Table

Item #	UOM
13000	1/cs

- The mast of the table pivots and can be locked in one of the three positions from flat on the floor to 90°. (Figure A)
- New and improved larger tubing.
- The pivot feature allows the table top to be positioned closer to individuals in a bed, sitting in a wheelchair, or folded for storage.
- The tabletop tilts with the angle of the mast to provide a flat surface regardless of the angle of the mast.
- Four 2" casters (two locking).
- Only wheels have to be installed.

### specifications

Tabletop Adjustment	28" - 37"
Tabletop Surface Dimensions	15.75" (D) x 23.5" (W)
Overall Dimensions at Base	25" (W) x 28" (D)
Weight	21 lbs.
Carton Shipping Weight	22 lbs.
Weight Limit	40 lbs.
Warranty	Limited Lifetime

## Seat Lift Chair Table

Item #	Product Description	UOM
13085RN	Right	1/cs
13085LN	Left	1/cs

- Designed for use with a lift chair, standard recliner or couch.
- Tabletop swivels to allow for easy access and departure. (Figure A)
- Can be used for many activities such as eating, writing, crafts, working on the computer and reading.
- Right or left positioning with 180° rotation accommodates seating near any wall.
- Tabletop pivots 360° at the touch of a hand.
- Teak wood grain laminate finish is attractive.
- A built-in magazine rack provides added convenience.

### specifications

Tabletop Adjustment	23.5" - 33"
Tabletop Surface Dimensions	15" (D) x 30" (W)
Overall Dimensions at Base	N/A
Weight	41 lbs.
Carton Shipping Weight	44 lbs.
Weight Limit	50 lbs.
Warranty	One Year Limited