## patient room



















## **Economy IV Pole** 4-Leg with Removable Top

Item #	Product Size	Color	UOM	
13029	4 Hook		1/cs	
13033	2 Hook		1/cs	
13033SV	2 Hook	Silver Vein	1/cs	

13033

- 4 rubber wheels provide for a smooth transport and maneuverability.
- Removable hooks make this economy IV pole versatile.
- Easily converted to 2-hook or 4-hook with easy to release pushpin.
- Chrome plated steel with weighted base provides strength, durability and reduces risk of tipping.
- Height range of 40"-82".
- Locking collar allows for easy height adjustment.

## IV Pole Top (Only)

Item #	Product Size	UOM
13028	4 Hook	1/cs
13040	2 Hook	1/cs

- Easily attaches to 13016, 13029, 13033 and 13043 IV
- Chrome-plated steel is easy to maintain.

specifications	
Warranty	Limited Lifetime

specifications	
Warranty	Limited Lifetime