



















Premium Series Shower Chair with Back and Arms

Item #	Product Description	UOM
RTI 12505	With Back and Arms	4/cs

- Easy, tool-free height adjustment with height indicator. Adjusts in ½" increments. (Figure A)
- Enhanced comfort seat fits all shapes and sizes.
- Easy, tool-free removable back, legs and arms.
- Lightweight design has large stable legs with soft non-slip feet.
- Retail packaged.

Shower Chair with Folding Back

Item #	UOM
RTL12606	1/ea

- Adds safety and support in the bath tub and shower.
- Easy, tool-free assembly and height adjustable legs.
- Folding back (Figure A) and removable legs make this shower chair great for storage and transport.
- · Lightweight design with stable legs and non-slip feet.
- Add comfort and convenience with optional arms and soap trays. (Figure B)
- · Retail packaged.

specifications	
Outside Legs	17" (W) x 15.5" (D)
Seat Dimensions	21.75" (W) x 16.5" (D)
Seat Height	16" (L) × 20.5" (H)
Weight	10 lbs.
Carton Shipping Weight	48 lbs.
Weight Limit	350 lbs.
Warranty	Limited lifetime

specifications	
Outside Legs	17.375" — 18.875"
Seat Dimensions	21.5" (W) x 13.125"(L)
Seat Height	17"-23"
Weight	8 lbs.
Carton Shipping Weight	8.6 lbs.
Weight Limit	300 lbs.
Warranty	Limited lifetime