



Shower Stool with Padded Rotating Seat

Item #	UOM
RTL12061	4 retail boxes/cs

- The rotating padded seat reduces twisting and reaching while in the shower.
- Seat rotates 360° and locks every 90°. (Figure A)
- The rotating top of the seat can be removed for easy cleaning.
- Easy tool-free assembly.
- Comes with removable tray for storing personal items.
- Ideal for individuals who have balance or mobility issues.
- Seat height can be adjusted from 16" to 23", in 1" increments.
- Retail packaged.

specifications	
Outside Legs	15" (W) x 15" (D)
Seat Dimensions	14" (W) x 14" (D)
Seat Height	16" - 23"
Weight	5 lbs.
Carton Shipping Weight	34 lbs.
Weight Limit	450 lbs.
Warranty	Limited Lifetime

Swivel Seat Shower Stool

Item #	UOM
RTL12061M	4 retail boxes/cs

- The rotating padded seat reduces twisting and reaching while in the shower.
- Seat rotates 360° and locks every 90°. (Figure A)
- The rotating top of the seat can be removed for easy cleaning.
- Easy tool-free assembly.
- Comes with removable tray for storing personal items.
- Ideal for individuals who have balance or mobility issues.
- Seat height can be adjusted from 16" to 22", in 1" increments.
- Retail packaged.

specifications	
Outside Legs	15" (W) x 15" (D)
Seat Dimensions	14" (W) x 14" (D)
Seat Height	16" - 22"
Weight	6.6 lbs.
Carton Shipping Weight	28.65 lbs.
Weight Limit	450 lbs.
Warranty	Limited Lifetime