

## Deluxe Aluminum Shower Bench without Back

HCPCS: E0245




FIGURE A




FIGURE B



FIGURE C

- Tool-free assembly of seat and legs (Figure A)
- Drainage holes in bench reduce slipping
- Frame is lightweight, durable and corrosion proof
- New pinch-free cover allows for push pins to be depressed without pinching finger The audible "click" will be heard when secure (Figure B)
- Legs are height adjustable in 1" increments, and angled with suction style tips to provide additional stability
- Support collar prevents leg movement and reduces rattle
- Made using recycled materials

PRODUCTS	Description
12203KD-4	Bath Seat w/o Back, KD, 4/cs
SPECIFICATIONS	Description
Carton Shipping Weight	7 lbs.   7 lbs.   21 lbs.
Outside Legs	19.5" (W) x 19.5" (D)
Seat Dimensions	19.25" (W) x 11.5" (D)
Seat Height	15.5" - 19.5"
Warranty	Limited Lifetime
Weight	5 lbs.   5 lbs.   4.5 lbs.
Weight Capacity	300 lbs
ACCESSORIES	Description
12002T	Tips for 12002,4/pkg., BX