

















Tub Rail

Item #	Product Size	UOM	
12034	8" x 11"	1/bx	

- Durable, chrome-plated steel is attractive and easy-toclean.
- Rubber pads on all tub contact points eliminate scratches.
- Knurling around grab bar diameter provides a no-slip grip.
- Easy to install.
- Tools included for installation.

specifications	
Adjustable Clamp	3.5" (smallest) 5.5" (widest)
Height	Not adjustable
Weight	4.5 lbs.
Carton Shipping Weight	6 lbs.
Weight Limit	300 lbs.
Warranty	Limited Lifetime

Clamp-On Tub Rail

Item #	Product Description	UOM
RTL12036-ADJ	Adjustable Height, Parallel to the Tub	1/bx
12036	Non-Adjustable Height, Parallel to the Tub	1/bx
12041	Non-Adjustable Height, Perpendicular to the Tub	1/bx

- Durable, white powder-coated steel construction.
- 14.5" 17" height adjustment (for RTL12036-ADJ only).
- Width adjusts from 3" 7".
- Stainless steel locking mechanism and non-wearing rubber pads provide scratch-proof security.
- Fits most tubs (not for use with fiberglass tubs).
- Tool-free installation.

specifications	
Adjustable Clamp	3" (smallest) 7" (widest)
Height (RTL12036-ADJ)	14.5"- 17"
Height (12036 12041)	14.5"
Weight (RTL12036-ADJ)	7.0 lbs.
Weight (12036 12041)	6.5 lbs.
Carton Shipping Weight (12036 12041)	8.5 lbs.
Carton Shipping Weight (RTL12036-ADJ)	9 lbs.
Weight Limit	300 lbs.
Warranty	Limited Lifetime