



JACUZZI
RECLINING BACK SUPPORT
FOR WIDER TUBS

Model no. 1574456

Aquatec J

Model no. 1574455 AQUATEC J, Reclining Back WIDE - White
Model no. 1574456 AQUATEC J, Reclining Back WIDE - Blue

The Aquatec J is designed with extra wide side flaps to fit into Jacuzzi-style tubs. With the ability to recline to 40 degrees, this model provides the safety, comfort and ease of use the Aquatec series is known for.

Features

- Extra wide sideflaps to fit into Jacuzzi-style tubs
- Stable, robust frame
- Can be used with or without cover materials
- Ergonomic floating hand control
- Lithium-ion rechargeable battery pack
- Lift function will not allow the lift to descend if there is not enough battery power present to raise the user
- Can be easily folded and disassembled into two parts with one hand
- Self-releasing suction cups
- 40 degrees of recline

Specifications

Seat Depth	19"
Seat Width	15"
Seat Width with Side Flaps Down	34"
Overall Length	24", 33" reclined
Back Height	26"
Back Width	14"
Base Width	11"
Base Length	23"
Minimum Seat Height	2 ½"
Maximum Lift Height (Optional Risers up to 19 ¾")	17"
Weight Capacity (lbs.)	300
Lift Cycles per Charge	6