



R *ECLINE*
BACK SUPPORT

Model no. 1573972

Aquatec *R*

Model no. 1573972
Model no. 1573877

AQUATEC R, Reclining Bath Lift - Blue
AQUATEC R, Reclining Bath Lift - White

The Aquatec *R* reclining bath lift combines the ability to recline to 40 degrees with industry-leading 2 ½" minimum seat height. The floating hand pendant with spiral cable ensures that the control is always within reach when needed.

Features

- Stable, robust frame that separates easily into two parts
- Reinforced lifting mechanism
- Can be used with or without cover materials
- Ergonomic floating hand control
- Lithium-ion rechargeable battery pack
- Lift function will not allow the lift to descend if there is not enough battery power present to raise the user
- Self-releasing suction cups
- Can be easily folded and disassembled into two parts with one hand
- High back with integrated neck support
- 40 degrees of recline

Specifications

Seat Depth	19"
Seat Width	15"
Seat Width with Side Flaps Down	28"
Overall Length	24", 33" reclined
Back Height	26"
Back Width	14"
Base Width	11"
Base Length	23"
Minimum Seat Height	2 ½"
Maximum Lift Height (Optional Risers up to 19 ¾")	17"
Weight Capacity (lbs.)	300
Lifts Cycles per Charge	6