# commodes





orve

## Deluxe Bariatric Commode Assembled

Item #	Color	UOM
11130-2	Gray	2/cs

- 650 lb. weight capacity.
- Large, padded backrest provides patient comfort.
- Durable, heavy-duty, gray powder-coated steel tubing is sturdy and easy to maintain.
- Large, durable snap-on seat.
- Plastic armrests provide comfort and stability.
- Pail easily slides in from the front of the commode.
  Comes complete with 12 qt. commode bucket and
- cover that removes from behind.

### specifications

Width		23" between arms   27" Outside
Seat Dimensions		23.25" (W) × 18.5" (D)
Height		15.5" - 22"
Weight		17 lbs.
Carton Shipping Weight		42 lbs.
Weight Limit		650 lbs.
Warranty		Limited Lifetime
accessories		
RTL12085	Biodegradable Sanitary mode Liner	Com- 12 liners/box, 6bx/cs

## Deluxe Bariatric Drop-Arm Commode Assembled

ltem #	UOM
11135-1	1/cs

- 1000 lb. weight capacity.
- Easy-to-release arm mechanism allows for safe lateral patient transfers to and from commode.
- Durable, heavy-duty, gray powder-coated steel tubing is sturdy and easy to maintain.
- Large, durable snap-on seat.
- Padded armrests provide comfort and stability.
- Pail easily slides in from the front of the commode.
- Comes complete with 12 qt. commode bucket and cover that removes from front.

#### specifications

Width	25.75" betwee	n arms   26.25" Outside
Seat Dimension	S	23.25" (W) x 18.5" (D)
Height		17.5" - 22"
Weight		21 lbs.
Carton Shipping	g Weight	26 lbs.
Weight Limit		1000 lbs.
Warranty		Limited Lifetime
accessories		
RTL12085	Biodegradable Sanitary Com- mode Liner	- 12 liners/box, 6bx/cs