

300











Raised Toilet Seat with/without Lid

Item #	Product Description	Product Weight	Product Carton Shipping Weight	UOM
12062	2", Without Lid	2 lbs.	4 lbs.	1/cs
12063	2", With Lid	2.5 lbs.	4.5 lbs.	1/cs
RTL12064	4", Without Lid	3 lbs.	5 lbs.	1/cs
12065	4", With Lid	3.5 lbs.	5.5 lbs.	1/cs
12066	6", Without Lid	5 lbs.	6.5 lbs.	1/cs
12067	6". With Lid	5.5 lbs.	7 lbs.	1/cs

- Designed for individuals who have difficulty sitting down or standing up from the toilet.
- · Heavy-duty molded plastic construction.
- Locking device with larger, heavy-duty "worm screw" and locking plate for a safe, secure fit. (Figure A)
- Lightweight and portable. Fits most toilets.
- · No tools required for installation.
- · Easy to clean.

specifications	
Seat Dimensions	14" (W) × 16" (D)
Weight Limit	400 lbs.
Warranty	Limited Lifetime

Open, Padded Raised Toilet Seat with Four Locking Brackets

Item #	UOM
12030	1/cs

- Designed for individuals who have difficulty sitting down or standing up from the toilet.
- Adds 5" to height of toilet seat.
- Padded seat provides extra comfort.
- Open seat promotes patient independence.
- Four adjustable locking clamps provide versatility.
- Adjustable angle and seat height.
- Lightweight and portable.
- Easily attaches to toilet bowl and locks in place with two rear locks.
- · No tools required for installation.
- Easy to clean.
- · Fits most toilets.

specifications	
Seat Dimensions	16" (W) x 15.5" (D) x 5" (H)
Weight	6 lbs.
Carton Shipping Weight	9 lbs.
Weight Limit	300 lbs.
Warranty	Limited Lifetime