

Cefar® Basic



Cefar Basic

A pain relief stimulator meeting all the basic needs

Cefar Basic offers TENS programs for most of the common pain indications. Cefar Basic has 3 preset treatment programs: TENS Conventional, TENS Burst and Modulated TENS (pulse width) programs.

Features

- Cefar Easy Touch™: Automatic Intensity Adjustment - At the press of a button, the intensity is automatically set and adjusted according to the results of the Auto test analysis. Cefar Easy Touch™ makes treatment easier and ensures stimulation takes place at an optimal and effective level
- Backlit Screen: The Cefar Basic has a larger screen and the backlight makes it much easier to read, even if the surrounding light is poor. Other advantages: pause function, program lock, key lock, compliance, self check
- Two year warranty



Ordering Information

Part Number Description

Cefar Basic

1490-EXP	Cefar Basic, Complete kit (DE-EN-ES-FR-IT-NL-NO-SE)
----------	---

Standard Accessories

1240	Set of Lead Wires (2pcs)
1567	Cefar Basic User Manual (DE-EN-ES-FR-IT-NL-NO-SE)
1559	Cefar Bag (Black), 2ch
1495	Cefar Neck Strap
1551	Batteries, Alk 1.5V AA (2 pcs)

Optional Accessories

1570	Fast Charger 1.2V
1575	Rechargeable Batteries 1800 mAh

Technical Specifications

Number of Channels:	2
Number of Programs:	3 preset
Constant Current:	Up to a resistance of Max 1500 Ohm/channel
Maximum Amplitude:	60 mA
Intensity Increments:	Minimum 0.5 mA
Max Pulse Duration/Width:	180 µ seconds
Frequency:	80 Hz
Power Source:	2x1.5V AA
Dimensions*:	120 x 30 x 50 mm
Weight:	180g, including battery
Electrical Safety Class:	Class II, Type BF
Safety Tests:	IEC 60601-1, IEC 60101-1-2

* Dimensions are listed as L x W x H

Waveforms



- Asymmetrical Biphasic 100% compensated