

## Airgo eXcursion X23 Lightweight Side-fold Rollator

HCPCS:



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E

- Ultra-lightweight side-fold frame designed for superior stability
- Patented ergonomic handgrips absorb vibrations, while also placing your hand at a natural angle for better brake access
- Easy to fold by just a pull of the strap, and the covered soft basket folds with the frame
- Includes a folding clip for safer lifting and transporting
- Padded seat and full contoured backrest for superior comfort
- Locks in open position for added safety
- Reflective trim for added visibility in the dark
- Easy to activate hand brakes with included parking brake function
- 8" (20 cm) comfort-grip wheels with low-profile treaded tires positioned outside the handles for increased stability
- Tool-free brake adjustment
- Reflective soft basket can be removed
- BONUS curb climber included
- AADL# C225
- ADP# MW3GM0065
- Recognized by CAOT

### PRODUCTS

### Description

700-917 Airgo eXcursion Tall 23", P. Blue, 1 EA

700-917BK Airgo eXcursion Tall 23", Black, 1 EA

### SPECIFICATIONS

### Description

Handle Height 33.5" to 39.5"

Seat Height 23"

Seat Width 17.5"

Unit Width 25"

Warranty Lifetime on Frame, 5 Years on Handle & Brakes, All other Parts 3 Months

Weight 17.1 lbs.

Weight Capacity 300 lbs

### ACCESSORIES

### Description

730-942 Slow brakes sfold, Black, 1PR, 25PR/CV, PR

730-946 CANE HOLDER SFOLD, BLACK, 1EA, 100CV, EA

730-944 CUP HOLDER SIDE FOLD BLACK, 1EA, 50CV, EA

730-950 Apollo LED Light, 1EA, 48CV, Ret, EA