



## Combination Padded Transfer Bench/Commode Tool-free Removable Back and Arm



FIGURE A



FIGURE B



FIGURE C



FIGURE D

- Combines transfer bench and commode into one
- Tool-free assembly of back, legs and arm (Figure A - C)
- Comfortable cushioned seat and backrest
- Aluminum A-frame construction is lightweight and provides additional stability
- Height adjusts in ½" increments with unique “Dual Column” extension legs
- No exposed hardware to injure patient
- Extra-large suction cups provide added safety
- New pinch-free cover allows for push pins to be depressed without pinching finger (Figure D)

PRODUCTS	Description
12005KDC-1	Combination Padded Transfer Bench/Commode, 1/ea

  

SPECIFICATIONS	Description
Outside Legs	17.25" (W) x 28.5" (D)
Seat Dimensions	24" (W) x 16" (D)
Seat Height	18.5" - 22.5"
Weight	10 lbs.
Carton Shipping Weight	15 lbs.
Weight Limit	400 lbs.
Warranty	Limited Lifetime

  

ACCESSORIES	Description
12011SS	Sm Suction Cups Pr., 1pr/bx
12011SL	Lg Suction Cups Pr., 1pr/bx