

Vacu-Aide® Suction Unit

HCPCS: E0600



FIGURE A



FIGURE B

- The DeVilbiss 7305 Series Vacu-Aide® Suction Unit meets suction and capacity requirements in a small, lightweight design
- High Performance Vacuum adjustments allow for 80-550 mm Hg and a free flow of 27 LPM
- Adjustable Flow Regulator Easy-to-turn knob conveniently placed directly below gauge for easy use and control Knob has safety lock feature which will not allow the knob to be twisted off (Figure A)
- The Vacu-Aide® Suction Unit 7305P series meets the RTCA/DO160D Aircraft Standard as well as the ISO-10079 1:1999 Suction Standard
- Optional Battery Capabilities For use up to one hour Provides the patient with added security when needed most
- The high flow capabilities (27 LPM) of the Vacu-Aide® Suction Unit can return a patient to a state of comfort quickly
- Unit can be operated while inside the carrying case (Figure B)

PRODUCTS	Description
7305D-D	Homecare Suction Unit, NAM
7305D-D-EXF	Homecare Suction Unit, NAM 800cc EXF
7305P-D	Homecare Suction Unit w/ Btty, NAM
7305P-D-EXF	Homecare Sucti.Unit w/Btty,NAM 800cc EXF
SPECIFICATIONS	Description
Dimensions	9.0" (H) x 7.0" (W) x 8.0" (D)
Electrical Requirements	100-240 V AC 50/60 Hz .75 A max; 12 V DC, 33 W max
Flow Range	27 LPM Free Flow, Unrestricted
Vacuum Range	80-550 mm Hg
Warranty	Limited 2 Years on Unit, Limited 90 Days on Battery
Weight	6.3 lbs with Battery; 3.8 lbs without Battery
ACCESSORIES	Description
7305D-632	Int Fltr Disp Canstr PK48, CS
7305D-633	Int Filter Disp Canister Comp Kit Assmd, EA
7305D-633-25	Int Fltr Disp Cnstr Comp Kit Assmd Cs 25, CS
PARTS	Description
DV51D-606	US POWER CORD, EA
7304D-619	BLACK 12V DC PWR COR, EA
7314P-613	AC TO DC ADAPTER SEP SALES PKG, EA
099HD-614	LINE CORD PKG (dft), EA
7305P-631	CONTINENTAL EUROPE CORD, EA
4650D-609	PUSH-IN FEET (4-PK), PKG

PARTS	Description
7304D-636	VACUUM GAUGE TEST FIXTURE, EA
7305D-605	SEP-SALES COMPRESSOR ASS'Y, EA
7305P-614	SEP SALES BATTERY ASSY,1/ea, EA
7305D-614	ElbowConct,Filter, 3.375" tubing, 1ea/ea, EA
7305D-608	BACTERIA FILTER (NON-STERILE) 12EA/PKG, PKG
7305D-612	REPLACEMENT TUBING 6-PK, PKG
7305D-633	Int Filter Disp Canister Comp Kit Assmd, EA
7305D-633-25	Int Fltr Disp Cnstr Comp Kit Assmd Cs 25, CS
7305D-635	800ML DISP CONTR FILTER Cart. Asml,12 P, CS
7305D-619	ELBOW (DFT) 6-PK SS, PKG
7305D-632	Int Fltr Disp Canstr PK48, CS