



ROHO® Planar Solid Seat Insert

The ROHO Planar Solid Seat Insert provides a solid base for use with ROHO cushioning products. The base provides support to assist with positioning. Made of polyurethane coated birch plywood, the 3/8 in. (.3 cm) base easily fits under the cushion, inside the ROHO cushion cover.

Eliminates the hammock effect of sling seats

Limited warranty: 12 months



ROHO® Contour Base

The lightweight, molded base provides a solid foundation for sling seat wheelchairs. Easily fits under the cushion inside the ROHO® cushion cover. The base provides support to assist with positioning.

Limited warranty: 12 months



ROHO® Privacy Shield*

The ROHO Privacy Shield is perfect for anyone who wants coverage between the bottom of the back of their wheelchair and the seat base of their wheelchair.

Unlike other privacy shields on the market, the ROHO Privacy Shield is designed to attach directly to the canes, so it will work with most wheelchairs. It is easy to attach and detach for cleaning.

Limited warranty: 6 months

*This is a non-medical product.



Hand Pump, Single or Dual Action

ROHO Hand Inflation Pumps are compatible with all air-filled ROHO medical and LTV® products.