

### Features

High impact bed end panels are more durable, impact and scratch resistant; washable for easy cleaning

Universal bed ends can be used as either the headboard or footboard on the bed

Rubber caster tread for better tracking with non-marring material

DC low voltage, compact motor system  
Quiet motor operation and fast cycle time on lowering height of bed for efficient bed set-up


Motor system with waterproof enclosures is rated IP66. Washable with standard water pressure. The bed must not be power washed

Color-coded connections on junction box for easier maintenance

Double-insulated junction box provides extra protection against electrical shock

Ergonomic hand pendant comes with durable strain relief and heavy-duty cord. Waterproof casing for easy cleaning, conveniently attaches to bed sheet with clip

Interchangeability simplifies storage, reduces inventory complexity



1. Patented universal bed ends
2. Ergonomic hand pendant
3. Color coded connections
4. Emergency manual crank
5. Thermal overload protection
6. Double insulated electrical components/wiring

Model #	Dimensions	Sleep Surface	Weight Capacity	Shipping Weight	Carton Dimensions	Regulation Compliance	Limited Warranty
541 IVC	88" L 15"-23" H 36" W	80" L 36" W	450 lb.			UL962 CSA Certified	5 Years Welds 2 Years Mechanical/ Electrical
5000IVC Headspring				45 lb.	36" L 44" W 8" D		
549 IVC Footspring				81 lb.	37" L 47" W 8" D		
5301IVC Universal Bed End				55 lb.	25" L 39" W		