









## **PrimeMat 2.0 Impact Reduction Fall Mats**

Item #	Product Description	UOM
PM20BA	Fall Mat 2.0 - Brown	1/ea, 5/cs
PM20GA	Fall Mat 2.0 - Gray	1/ea, 5/cs

- Full 72" length provides maximum coverage alongside bed.
- Tough and lightweight Mats are constructed from high quality molded Polyurethane foam inner core (210kgs/cubic meter) to provide maximum absorption and rebound to reduce impact in the event of a fall.
- Passed 5 foot "Egg Drop Test" an egg was dropped from 5 feet above the mat without damage.
- Low profile with tapered edges allow equipment such as overbed tables or patient lifts to roll over easily.
- Anti-slip bottom netting keeps the mat secure to help prevent accidental tripping.
- Embossed non-slip PVC leather look top is not only attractive but extremely durable and puncture resistant.
- Non-slip pattern helps minimize slipping and is easy to clean with mild soap and water.
- Available in two colors; Brown & Gray.
- 1 Year Warranty.

	•	ere .	• •	
sn	eci.	tica:	tion	•
-	-	II Cu	tions	۱

Unit Dimensions 72" x 24" x 75"