

# Fall Sensor Pads With No Cords



**TL-2100G**  
Ghost Cord®  
Fall Monitor

**No Cord!**  
*The transmitter is in the pad!*



**TM-02**  
transmitter inside  
the neck of the pad.  
**PATENT PENDING**

## Go Cordless

The Smart Ghost Cord® system frees caregivers from managing cords between sensor pads and the alarm monitor allowing more options for monitor placement\* (even outside the room) as well as reducing tripping hazards.

Moving the alarming monitor outside of the room has two primary benefits:

- Staff members hear the alarm better
- Residents are not disturbed by alarm noise

MFR. #	SMART PRODUCT DESCRIPTION
TL-2100G	Ghost Cord® Wireless Pad Alert Monitor (1 year warranty)
TM-02	Ghost Cord® Wireless Pad Alert Transmitter Insert
GB-WI	Ghost Cord® Pad for Bed - 20"x 30" (1 year warranty)
GB-RI	Ghost Cord® Pad for Bed - 10"x 30" (1 year warranty)
GB-45	Ghost Cord® Pad for Bed - 10"x 30" (45 day warranty)
GB-90	Ghost Cord® Pad for Bed - 10"x 30" (90 day warranty)
GC-WI	Ghost Cord® Pad for Chair - 10"x 15" (1 year warranty)
GC-45	Ghost Cord® Pad for Chair - 10"x 15" (45 day warranty)
GC-90	Ghost Cord® Pad for Chair - 10"x 15" (90 day warranty)
FM-07C	Cordless Extra Long Floor Mat - 24"x 48" (1 year warranty)

\*up to 50 feet



**FM-07C Cordless**  
Floor Mat and  
**GC-45 Ghost Cord®**  
Sensor Pad  
**PATENT PENDING**

## Get Paged immediately from Cordless Pads

Caregivers can now be paged directly from a cordless sensor pad to a pager carried in a pocket when using a 433-CMU central monitor – no alarm monitor needed! Perfect for a quiet facility.

**433-PGD**  
Pager w/ LCD Display



**New!**  
**Paging**



**433-CMU**  
Central Monitoring Unit

**SMART CAREGIVER CORPORATION**

www.westcoastmedicalsupply.com 1-866-517-8243

SS-TL-2100G-REV.D

