



The Invacare I-Class Paddle Walker features a wide, deep frame for user convenience and is lightweight, making it easy to lift and maneuver. The paddle release mechanism makes folding the walker easy and is ideal for users with dexterity issues.

Composite lower side brace for added stability  
Folding mechanisms are easy to reach from a seated position  
Lighter and easier to lift and maneuver  
Large range of height adjustments  
PVC handgrips wear longer and provide greater comfort

**Warranty:** Limited Warranty: 5 Year

**Width Inside Hand Grips:** 17"

**Width Inside Base Legs:** 20.5"

**Overall Depth:** Opened 16.5" Folded: 4"

**Height Adjustment:** 25.4" - 32.4"

**Patient Height:** 4'4" - 5'7"

**Product Weight:** 5 lb. 1 oz.

**Product Weight Capacity:** 300 lb.