

Item Number	Description	Wheels	Weight Capacity Lbs (kg)	Approximate User Height Inch (cm)	Width Inside Hand Grip Inch (cm)	Height Adjustment Inch (cm)	Folded Width Inch (cm)	Pkg
Premium Two-Button Walker								
G07755	Adult	None	350 (159)	5'5"-6'4" (165-193)	16.5 (42)	32.5-39.5 (82-100)	4.5 (11)	4/cs
G07756	Junior	None	350 (159)	5'-6' (152-183)	16.5 (42)	27.25-34.25 (69-87)	4.5 (11)	4/cs
G30755P	Adult	None	300 (136)	5'5"-6'4" (165-193)	16.5 (42)	32.5-39.5 (82-100)	4 (10)	4/cs
G30756P	Junior	None	300 (136)	5'-6' (152-183)	16.5 (42)	27.5-34.5 (70-88)	4 (10)	4/cs
G30757W	Adult	5 (13)	300 (136)	5'5"-6'4" (165-193)	16.5 (42)	32.5-39.5 (82-100)	4 (10)	4/cs
G30758W	Junior	5 (13)	300 (136)	5'-6' (152-183)	16.5 (42)	27.5-34.5 (70-88)	4 (10)	4/cs
Standard Two-Button Walkers								
MDS864104	Adult	None	300 (136)	5'5"-6'4" (165-193)	16.75 (43)	32-39 (81-99)	4.5 (11)	4/cs
MDS86410W54	Adult	5 (13)	300 (136)	5'5"-6'4" (165-193)	16.75 (43)	32.5-39.5 (82-100)	4.5 (11)	4/cs
MDS86410W4	Adult	3 (8)	300 (136)	5'5"-6'4" (165-193)	16.75 (43)	32-39 (81-99)	4.5 (11)	4/cs
MDS86410J4	Junior	None	300 (136)	4'6"-5'5" (137-165)	16.75 (43)	25-32 (64-81)	4.5 (11)	4/cs
MDS86410JW54	Junior	5 (13)	300 (136)	4'6"-5'5" (137-165)	16.75 (43)	25.5-32.5 (65-83)	4.5 (11)	4/cs
MDS86410WJ4	Junior	3 (8)	300 (136)	4'6"-5'5" (137-165)	16.75 (43)	25-32 (64-81)	4.5 (11)	4/cs
Basic Two-Button Walker								
MDS864104B	Adult	None	400 (181)	5'5"-6'4" (165-193)	16.75 (43)	32-39 (81-99)	4.5 (11)	4/cs
MDS86410W54B	Adult	5 (13)	400 (181)	5'5"-6'4" (165-193)	16.75 (43)	32.5-39.5 (82-100)	4.5 (11)	4/cs
MDS86410W4B	Adult	3 (8)	400 (181)	5'5"-6'4" (165-193)	16.75 (43)	32-39 (81-99)	4.5 (11)	4/cs
MDS86410J4B	Junior	None	400 (181)	4'6"-5'5" (137-165)	16.75 (43)	25-32 (64-81)	4.5 (11)	4/cs
MDS86410JW54B	Junior	5 (13)	400 (136)	4'6"-5'5" (137-165)	16.75 (43)	25.5-32.5 (65-82.5)	4.5 (11)	4/cs
MDS86410JW4B	Junior	3 (8)	400 (136)	5'5"-6'4" (165-193)	16.75 (43)	25-32 (64-81)	4.5 (11)	4/cs

Item Number	Description	Wheels	Weight Cap Lbs (kg)	Approximate User Height Inch (cm)	Width Inside Hand Grip Inch (cm)	Height Adjustment Inch (cm)	Folded Width Inch (cm)	Pkg
Premium One-Button Folding Walker								
G30760P	Adult	None	300 (136)	5'5"-6'4" (165-193)	17.25 (44)	31.5-38.5 (80-98)	4.5 (11)	4/cs
G30761P	Junior	None	300 (136)	5'-6' (152-183)	17.25 (44)	27.25-34.25 (69-87)	4.5 (11)	4/cs
Standard One-Button Folding Walker								
MDS86615	Adult	None	250 (113)	5'-6' (152-183)	18 (46)	32-36 (81-91)	4.5 (11)	2/cs
MDS86617	Junior	None	250 (113)	4'6"-5'5" (137-165)	18 (46)	28-32 (71-81)	4.5 (11)	2/cs