



# SEAT CANES

Folds for Easy Transportation and Storage • Supports from 200 to 250 lbs.

**9-0227-00**  
Aluminum four leg seat cane with canvas saddle seat. Seat height is 21" and supports up to 250 lbs.




**9-0228-00**  
Steel four leg black seat cane with 12" round padded seat. Seat height is 17-1/2" and supports up to 250 lbs.



**9-0224-00**  
Aluminum seat cane silver with three leg and square black seat. Seat height is 21" and supports up to 250 lbs.



**9-0225-00**  
Adjustable, aluminum three leg seat cane, bronze with black oval seat and supports up to 200 lbs.



**9-0229-00**  
Aluminum three leg seat cane bronze with black formed seat. Seat height is 20" and supports up to 200 lbs.