













Shower Chair with Back and Removable Padded Arms

| Item # | UOM |
|---------|------|
| 12445-2 | 2/cs |
| 12445-1 | 1/cs |

- Removable padded armrests provide comfort and added security.
- Blow-molded bench and back provide comfort and strength.
- Drainage holes in bench reduce slipping.
- Aluminum frame is lightweight, durable and corrosion proof.
- · Leg height adjusts in 1" increments.

| specifications | |
|----------------------------------|-----------------------|
| Outside Legs | 17.75 (W) x 19.25 (D) |
| Seat Width | 16" |
| Seat Depth | 14.5", 16.25" |
| Seat Height | 16" - 21" |
| Weight (-1 -2) | 9.75 lbs. 19 lbs. |
| Carton Shipping Weight (-1 -2) | 10 lbs. 21 lbs. |
| Width Between Arms | 18.5" |
| Weight Limit | 300 lbs. |
| Warranty | Limited Lifetime |
| | |

Folding Shower Chair

| Item # | Product Description | UOM |
|--------|----------------------------|------|
| 12486 | Without Back | 1/cs |
| 12487 | With Back | 1/cs |

- Product available with back (main image) and without. (Figure A)
- Conveniently folds flat. (Figure B)
- Blow-molded bench provides comfort and strength.
- Drainage holes in seat reduce slipping.
- Aluminum frame is lightweight, durable and corrosion proof.
- · Angled legs provide additional stability.
- Retail packaged.

| specifications | |
|--|---------------------------------|
| Outside Legs (12486) | 20.5" (W) x 19.25" (D) |
| Outside Legs (12487) | 17.75" (W) x 14" (D) |
| Seat Dimensions (12486) | 19.75" (W) x 12" (D) |
| Seat Dimensions (12487) | 18" (W) x 11" (D) |
| Seat Height (12486 12487) | 17.5" - 18.5" 16.25" - 18.25" |
| Weight (12486 12487) | 4.5 lbs. 5 lbs. |
| Carton Shipping Weight (12486 12487) | 6.6 lbs. 7 lbs. |
| Weight Limit | 300 lbs. |
| Warranty | Limited lifetime |