



Transfer Tub Bench

Knocked Down Tool-free Back, legs and Arms

Item #	Product Description	UOM
12011KD-2	Knocked Down, Tool-free Back, Legs, and Arms	2/cs
12011KD-1	Knocked Down, Tool-free Back, Legs, and Arms	1/cs
RTL12011KDR	Knocked Down, Tool-free Back, Legs, and Arms	1/cs

- Tool-free assembly of back, legs and arm. (Figure A - C)
- A-frame construction provides stability.
- Durable blow-molded plastic bench and backrest.
- Height adjusts in ½" increments with unique "Dual Column" extension legs.
- Reversible to accommodate any bathroom.
- Pinch-free lever allows for push pins to be depressed without pinching fingers.

specifications

Outside Legs	32" (W) x 23" (D)
Seat Dimensions	26" (W) x 18.5" - 19.5" (D)
Seat Height	17.5" - 22.5"
Weight	10 lbs.
Carton Shipping Weight (KDR -1 -2)	13 lbs. 13 lbs. 24 lbs.
Weight Limit	400 lbs.
Warranty	Limited Lifetime

Padded Transfer Bench

Tool-free Assembly Back, Legs and Arms

Item #	Product Description	UOM
12005KD-2	Knocked Down, Tool-free Assembly Back, Legs and Arms	2/cs
12005KD-1	Knocked Down, Tool-free Assembly Back, Legs and Arms	1/cs
12005KDR-1	Knocked Down, Tool-free Assembly Back, Legs and Arms	1/cs

- Comfortable cushioned seat and backrest.
- Tool-free assembly of back, legs and arm. (Figure A - C)
- New pinch-free cover allows for push pins to be depressed without pinching fingers.
- A-frame construction provides stability.
- Height adjusts in ½" increments with unique "Dual Column" extension legs.

specifications

Outside Legs	29" (W) x 19.75" (D)
Seat Dimensions	24" (W) x 16" (D)
Seat Height	17.75" (L) x 21.75" (H)
Weight	11.75 lbs.
Carton Shipping Weight (-1 -2)	14 lbs. 26.5 lbs.
Weight Limit	400 lbs.
Warranty	Limited Lifetime