

Padded Transfer Bench

HCPCS: E0247



FIGURE A



FIGURE B




FIGURE C

- Comfortable cushioned seat and backrest
- Tool-free assembly of back, legs and arm (Figure A - C)
- New pinch-free cover allows for push pins to be depressed without pinching fingers
- A-frame construction provides stability
- Height adjusts in ½" increments with unique "Dual Column" extension legs

PRODUCTS	Description
12005KD-1	Transfer Bench, Padded, KD, 1/cs
SPECIFICATIONS	Description
Carton Shipping Weight	14 lbs. 26.5 lbs.
Outside Legs	29" (W) x 19.75" (D)
Seat Dimensions	24" (W) x 16" (D)
Seat Height	17.75" (L) x 21.75" (H)
Warranty	Limited Lifetime
Weight	11.75 lbs.
Weight Capacity	400 lbs
ACCESSORIES	Description
12011SL	SuctionCups,Lg.12011,1pr/bx, BX
12011SS	SuctionCups,Sm.12011,1pr/bx, BX
PARTS	Description
12005KD-B	Back Assembly With Pad,12005KD,1ea, EA
12005KD-ARM	Arm,12005KD,1ea, EA
12005KDSP	Seat Pad,12005KD,1ea, EA
12005KD-LEGS	Legs With Tips,12005KD,Set of 4,1set/pkg, PKG
12011SL	SuctionCups,Lg.12011,1pr/bx, BX
12011SS	SuctionCups,Sm.12011,1pr/bx, BX
12005KDTP	Transfer Pad,12005KD and KDC,1/ea, EA