

## Accessories



### Anti-Tippers

*For additional safety.*

Model no. 1470833

Available for Ocean VIP only



### Reclining Backrest

*Two stainless steel gas springs make the step-less adjustment of the reclining backrest easy. Reclines 10°-50°.*

Model no. 1470085

Available for Ocean VIP only



### XL Backrest

*Increases the width of the backrest by 3".*

Model no. 1534979 (prior to SN 000014541064 Ocean XL/  
000014551780 Ocean)

Model no. 1558847 (after SN 000014541065 Ocean XL/  
000014551781 Ocean)

### XL Arms

*Adds an additional 3" between arms.*

Model no. 1470086

Available for Ocean VIP and Ocean E-VIP only




### Frontal Splash Guard

*Easily fitted. For soft seat with hygiene recess.*

Model no. 16511

Available for all Ocean models



### Armrest Lock

*Locks the armrest in the lower position.*

Model no. 16597

Available for Ocean Standard, SP, XL (up to SN 14403119) only

Model no. 1470077

Available for all Ocean Models (after SN 14403200)



### Chest Belt/Pelvis Belt

*Comfortable, soft, machine washable belts affix easily to enhance patient safety. Adjustable 26.5" - 41".*

Model no. 1470081, Chest Belt

Model no. 1470082, Pelvis Belt

Available for all Ocean models

### Belt Maintenance

**REGULAR CLEANING:** Machine wash, hot, 60°C (140°F) using laundry detergent. Tumble dry, or drip dry. Do not iron.

**SPECIAL CLEANING:** For additional sanitization, the belt can be machine washed, very hot, 70°C (160°F) with bleach.