



Mermaid Specifications	
Safe Working Load: 125 kg / 19.6 st	Maximum Seat Depth: 410 mm
Maximum Height to Base of Seat: 740 mm	Centre of Seat to Centre of Mast: 430 mm
Minimum Height to Base of Seat: 140 mm	Maximum Seat Turning Radius: 850 mm
Maximum Seat Travel: 485 mm	Maximum Clearance from Mast (to allow for seat rotation): 1000 mm
Maximum Seat Width: 450 mm	Total Weight - Mermaid Manual: 29 kg
Maximum Width Between Arms: 395 mm	Total Weight - Mermaid Electric: 37 kg