

ALBAC RESCUE MAT



Albac Rescue Mat is used all over the world for fast, safe evacuation of children or adults indoors and out, horizontal or vertical rescue. Glides easily over all types of terrain. Mat is heat resistant and ready for instant use. Stretcher like design allows single person operation, even on staircases. Thermal blanket enclosed for patient comfort.

NEW Bariatric Size Available

- ✓ Fluorescent straps for detection in dim lighting
- ✓ Foot rest pocket for comfort and security of patient
- ✓ Polypropylene board acts as insulator over surfaces, prevents friction burns to patients
- ✓ Eleven separate handles for rescuer's ease
- ✓ Simple instructions clearly visible on self-storage pocket

Advanced Egress Solutions, Inc.

Certified Woman Owned Business (WBE) New York State, New York City and Nassau County



Post Office Box 33
Bethpage, NY 11714
Phone: 516 359-8521

E-mail: sales@advancedegresssolutions.com
Web site address
www.advancedegresssolutions.com




The mat has been tested to hold over 1000 pounds. Flammability testing revealed slow to no ignition. Copies of all testing is included in our Owner's Manual. The bottom of the mat is made of a polypropylene board that slides easily, absorbing any concussions in transit. For patient comfort.

ALBSTD01 Standard Size / open 75"x 26", stored 14"x7"x7" ALBLRG02 Large Size / open 75"x 32", stored 14" x 7" x 7" ALBBAR09 Bariatric* Size / open 46" x 8' stored 25" x 9" x 7" ALBHLD03 AlbacMat Holder-neatly stores 4 Standard/Large Mats *Bariatric Size Mat requires a minimum of 2 handlers.

Albac has a foot rest pocket and thermal blanket included. The straps are fluorescent for easy detection in dim lighting. Great for both vertical and horizontal evacuation, indoors and out.

The AlbacMat Holder neatly stores 4 mats and can be carried to the site of an emergency for instant use.



The AlbacMat is used in Hospitals, Nursing Homes, Schools/Universities, Mining Industry, Military/Defense, Fire/Rescue/Emergency, Search and Rescue, Airline Industry, Worksite



Safety, Amusement Park Safety, Trains/Subway/Bus, High Rise Buildings, Stadiums/Arenas, Mortuaries, Cruise Ships, Private and Group Homes, Anywhere!

