

## Pulmo-Aide® Compressor Nebulizer System

HCPCS: E0570,A7003



- The Pulmo-Aide Compressor Nebulizer System is the only approved compressor for use with TOBI (Tobramycin)
- One of the longest lasting compressors on the market with an average product life of 10 years
- Unmatched quality and durability Great for nursing homes and long-term users
- Metal frame motor construction
- Large built-in compartment for storing nebulizers and tubing
- Built-in nebulizer holder is a convenient place to store your nebulizer if your treatment is interrupted

PRODUCTS	Description
5650C	PULMO-AIDE CANADA
SPECIFICATIONS	Description
Dimensions	6.3" (W) x 10.5" (H) x 10.1" (D)
Free Air Flow	12 lpm Typical
Max Pressure	30 psig or Greater
Sound Level	51-55 dBA
Warranty	Limited 5 Year
Weight	7.1 lbs
ACCESSORIES	Description
5650D-636	