















Item #	UOM
11101W-2	2/cs

- Easy-to-release arm mechanism allows for safe lateral patient transfers to and from commode.
- When arm is folded, it provides a handle for additional support while transferring. (Figure A - B)
- Rust-resistant, swivel casters (2 with locks).
- White powder-coated steel frame. Padded armrests provide extra comfort and support.
- Comes with 7.5 qt. commode bucket with carry handle and splash shield.

specifications	pecifications	
Width	18" between arms 22" Outside	
Seat Dimensions	14.25" (W) x 16.75" (D)	
Height	17" - 21"	
Weight	22 lbs.	
Carton Shipping Weight	47 lbs.	
Weight Limit	300 lbs.	
Warranty	Limited Lifetime	
3		

accessories		
RTL12085	Biodegradable Sanitary Commode Liner	12 liners/box, 6bx/cs

Aluminum Commode with Wheels

- Brass push button allows for easy, tool-free assembly and disassembly of back. (Figure A)
- Non-skid, rust-resistant, swivel casters (2 with locks). (Figure B)
- Aluminum construction is lightweight and durable.
- Angled back for support and comfort.
- Full length plastic arms for added comfort.
- Comes with a 7.5 qt. commode bucket, splash shield and cover. Splash shield acts as pail holder.

specifications	
Width	19.5" between arms 23.5" Outside
Seat Dimensions	13.75" (W) x 16.5" (D)
Height	19.25"
Weight	12 lbs.
Carton Shipping Weight	28 lbs.
Weight Limit	300 lbs.
Warranty	Limited Lifetime

accessor		
RTL12085	Biodegradable Sanitary Commode Liner	12 liners/box, 6bx/cs