

## Airgo Comfort-Plus Lightweight Transport Chair

HCPCS:



- New and improved with more security and comfort enhancing features
- Flip-back detachable arms make it easier to transfer in and out, or simply get closer to a table
- Comfort-Grip 8" wheels provide shock absorption and better brake performance
- Anatomically correct seat-to-armrest height for a naturally erect posture
- Fold-down back for easy transportation and storage
- Removable swing-away footrests with tool-free adjustment of length for greater comfort
- Heel-loops and seatbelt included as standard features for increased safety
- Dual locking rear wheels for safe, secure transfers
- Chip-resistant finishing and no sharp edges on hardware
- Lightweight and compact

PRODUCTS	Description
700-841	Airgo Comfort-Plus Lightweight Transport Chair, EA
700-846	Airgo Comfort-Plus Lightweight Transport Chair, EA
700-848	Airgo Comfort-Plus Lightweight Transport Chair, EA
SPECIFICATIONS	Description
Folded Width	10.5
Handle Height	37"
Seat Dimensions	17" or 19" (W) x 16" (D)
Seat Height	19"
Seat Width	17"   19"
Unit Dimensions	22" or 24" (W) x 26" (D)
Warranty	Frame and Cross-brace 3 Years, All Other Metal Parts 6 Months, Upholstery, Plastic and Rubber 3 Months
Weight	20.5 lbs. net
Weight Capacity	300 lbs