

ROHO® ADAPTOR® PAD

The ADAPTOR Pad turns any flat or contoured surface into a safe, cushioned, protective environment. Since the air cells are not linked to one another, the ADAPTOR can be custom-cut and fit to address a particular need.

It's not a seat; it's designed to protect other pressure points. It can be cut to fit arm rests, trays, foot rests or anywhere protection is needed.

- ✓ Individually sealed air cells
- ✓ Skin/soft tissue protection
- ✓ Can be cut to size
- ✓ Absorbs shock



Dimensions

Height: 1 in. / 2 cm
Width: 9.75 in. / 25 cm
Depth: 13.5 in. / 34 cm
Weight: 0.5 lbs. / 0.2 kg

Weight limit

None, if properly sized
Limited warranty
24 months

ROHO® HEAL PAD® Cushion

The HEAL PAD places your foot in an adjustable, therapeutic environment that reduces shear and friction on your heel. By supporting the most high-risk pressure points, the pad protects you from injury and facilitates healing.

- ✓ Adjustable
- ✓ Air-filled cells
- ✓ Available in singles and pairs



Dimensions

Height: 1.5 in. / 4 cm
Width: 11.5 in. / 29.5 cm
Depth: 12.5 in. / 31.5 cm
Weight: 0.75 lbs. / 0.9 kg

Limited warranty

24 months