

GUARDIAN[®] Bariatric Canes

- Designed with reinforced tubing and/or wider bases to achieve maximum weight capacity allowable
- High weight-bearing capabilities
- 500 lb (227 kg) Capacity



Bariatric Quad Canes

All steel construction with double-plated base allows for a 500 lb (227 kg) weight capacity. Contains a secure locking ring and dual-tubing insert for added strength and safety.

Single Point Bariatric Cane

All steel construction with dual-tubing insert for added strength and security. Single point design gives additional balance support and weight-bearing capabilities.

Item Number	Description	Color	Weight Cap Lbs (kg)	Adjustment Inch (cm)	Footprint Inch (cm)	Cane Weight Lbs (kg)	Pkg
Bariatric Quad Canes							
MDS86228XW	Large Base	Black	500 (227)	29-38 (74-96.5)	8 x 12 (20 x 30)	4.6 (2)	Each
MDS86222XW	Small Base	Black	500 (227)	29-38 (74-96.5)	6 x 8 (15 x 20)	3.6 (2)	Each
Single Point Bariatric Cane							
MDS86420XW	Single Point	Black	500 (227)	29-38 (74-96.5)	N/A	3.0 (1)	Each

Cane Tips and Grip

- Tips are reinforced with metal inserts for longer wear



Item Number	Description	Size: Inches (cm)	Color	Pkg
Cane Tips				
G30561-PC	Tip	5/8 (1.6)	Black	8 pr/cs
G30572-PC	Tip	3/4 (2.0)	Gray	8 pr/cs
G30581-PC	Tip	7/8 (2.0)	Black	8 pr/cs
G02162	Tip, Quad Cane	1/2 (1.3)	Gray	8 set/cs
G02172	Tip, Quad Cane	5/8 (1.6)	Gray	8 set/cs
MDS80212T	Tip, Replacement for MDS80212	7/8 (2.0)	Black	2 ea/pr
MDS80217T	Tip, Replacement for MDS80217	1 (2.5)	Black	2 ea/pr
MDS86424W	Tip, Quad Cane	1/2 (1.3)	Black	6 pr/cs
MDS86425W	Tip, Quad Cane	5/8 (1.6)	Black	6 pr/cs
MDS86426W	Tip, Quad Cane	3/4 (2.0)	Black	6 pr/cs
Cane Grip				
MDS86427W	Grip, for 420, 222, 228 Canes	4 (10.2)	Black	6 pr/cs