

## Hemi Sidesteppers and Non-Folding Walkers

- Lightweight design allows the hemi sidestepper to be maneuvered easily even with just one arm
- Hemi sidestepper and non-folding walker provide more stability and weight-bearing ability than a quad cane
- Hemi Sidestepper and non-folding walker are height-adjustable in 1" (2.5 cm) increments



G07770



MDS86620



MDS86204

### Premium Hemi Sidestepper

Lightweight frame and folding capabilities make this walking aid ideal for travel and storage. Bi-level hand grips offer extra assistance for users getting in and out of chairs.

### Standard Hemi Sidestepper

Wider base provides more stability than a quad cane. Sidestepper provides support and stability to those with limited or no dexterity in one arm or hand.

### Basic Non-Folding Walker

Secure frame without folding mechanisms eliminates even the slightest wobble for users needing maximum assistance with balance.

Item Number	Description	Wheels	Weight Cap Lbs (kg)	Width Inside Hand grip Inch (cm)	Height Adjustment Inch (cm)	Folded Width Inch (cm)	Pkg
<b>Premium Hemi Sidestepper</b>							
G07770	Premium Hemi, Adult	None	250 (113)	N/A	32-36 (81-91)	3 (7.6)	2/cs
G07771	Premium Hemi, Junior	None	250 (113)	N/A	28.5-32 (72-81)	3 (7.6)	2/cs
<b>Standard Hemi Sidestepper</b>							
MDS86620	Standard Hemi	None	250 (113)	N/A	30-35 (76-89)	3 (7.6)	Each
<b>Non-Folding Rigid Walker</b>							
MDS86204	Non-Folding Walker	None	250 (113)	18 (46)	31-37.5 (79-95)	N/A	4/cs