

## Wall Grab Bars

**GUARDIAN**

- Locate grab bars at any angle in the bathroom to help prevent falls
- Use grab bars in the kitchen, along stairways, or in the laundry room for additional fall prevention
- Combine grab bars with an elevated toilet seat or transfer bench for additional stability



### Premium Knurled Chrome

Adjustable oval flange rotates to allow for secure mounting. Knurled finish provides steady grip even when wet.



### Standard Knurled Chrome

Knurled-texture steel bars are easy to grip even in wet bathing areas.



### Basic White Finish

White painted steel with white screws blend in to match most bathroom decors.

Item Number	Length Inch (cm)	Diameter Inch (cm)	Pkg
G81020KE	12 (30.5)	1 (2.5)	3/cs
G81030KE	16 (40.6)	1 (2.5)	3/cs
G81040KE	18 (45.7)	1 (2.5)	3/cs
G81050KE	24 (61.0)	1 (2.5)	3/cs
G81060KE	32 (81.3)	1 (2.5)	3/cs
G81070KE	36 (91.4)	1 (2.5)	3/cs

Item Number	Length Inch (cm)	Diameter Inch (cm)	Pkg
MDS86012CHRG	12 (30.5)	1 (2.5)	3/cs
MDS86016CHRG	16 (40.6)	1 (2.5)	3/cs
MDS86018CHRG	18 (45.7)	1 (2.5)	3/cs
MDS86024CHRG	24 (61.0)	1 (2.5)	3/cs
MDS86032CHRG	32 (81.3)	1 (2.5)	3/cs

Item Number	Length Inch (cm)	Diameter Inch (cm)	Pkg
G87520	12 (30.5)	1 (2.5)	3/cs
G87530	16 (40.6)	1 (2.5)	3/cs
G87540	18 (45.7)	1 (2.5)	3/cs
G87550	24 (61.0)	1 (2.5)	3/cs
G87560	32 (81.3)	1 (2.5)	3/cs
G87570	36 (91.4)	1 (2.5)	3/cs

**Retail Packaging** Ready for Display

## Bath Tub Grab Bars and Safety Pole

- Give patients confidence and support with an extra grip to help them transfer into the tub



G98006

### Deluxe Tub Assist Bar

Step through clamp design provides unobstructed entry/exit path allowing better hand positioning. Tool-free installation has stainless steel components with rubber-lined pads that won't mar tub surface. Adjust to fit tubs 2.75" - 6.25" (7-16 cm).



G98080

### Bathtub Safety Rail

Clamps on two sides with foam lined clamp that won't mar tub surface. Mounting brackets adjust from 3.5" - 6" (9-15.24 cm) to fit most standard size tubs.



G98120P

Floor to ceiling pole that safely assists standing, sitting, climbing stairs, or transferring.

Item Number	Description	Pkg
G98006	Deluxe Tub Assist Bar	3/cs
G98080	Bathtub Safety Rail	3/cs
G98120B	Safe-T-Pole, Base	Each
G98120P	Safe-T-Pole, Pole	Each

**Retail Packaging** Ready for Display