

## Elevated Toilet Seats

- Makes using the toilet easier for those with trouble bending
- Available in a variety of locking mechanisms and mounting methods



MDS80326

### Extra Wide with Legs 600 lb (272 kg) Capacity

Steel legs go to the floor and soft foam armrests provide maximum support in getting on and off the toilet.



G30270A

### Locking with Padded Arms

Padded arm rests assist in getting on and off toilet. Locking device secures seat to toilet and requires no tools.



MDS80320S

### Locking with Molded Arms

Arms assist in getting on and off the seat. Locking device secures the seat to the toilet.



MDS80315

### Universal Locking

Height-adjustable brackets allow seat to fit over any toilet. Brackets are vinyl covered to prevent chipping or scratching of porcelain. Seat is height-adjustable with plastic splash guard to make cleaning easy.



MDS80324S

### Hinged Bolt-On








Seat can be lifted and lowered with hinged design. Existing seat and lid make this model less noticeable.



G30250

### Basic Raised Seat

Simple design made of polypropylene resin is easy to clean. Non-slip pads and lip inside toilet offer additional stability.

Item Number	Description	Wt Capacity Lbs (kg)	Seat Elevation Inch (cm)	Pkg
<b>With Arms</b>				
MDS80326	Extra Wide with Legs	600 (272)	4 (10)	Each 
G30270A	Locking with Padded Arms	350 (159)	5 (12)	3/cs 
MDS80320S	Locking with Molded Arms	300 (136)	5 (12)	Each 
<b>Without Arms</b>				
MDS80315	Universal Locking	250 (113)	4.75-6.75 (12-17)	Each
MDS80324S	Standard Hinged Bolt-On	300 (136)	4 (10)	Each 
MDS80324E	Elongated Hinged Bolt-On	300 (136)	4 (10)	Each 
G30260*	Locking Raised	350 (159)	5 (12)	3/cs 
G30250	Basic Raised Seat	300 (113)	5 (12)	3/cs 

\* Not Pictured

 **Retail Packaging** Ready for Display

## Toilet Safety Frame



G30300

Item Number	Wt Capacity Lbs (kg)	Height x Width Inch (cm)	Pkg
G30300	250 (113)	26-31 (66-78.7) x 18-24 (46-61)	2/cs
G30300-1	250 (113)	26-31 (66-78.7) x 18-24 (46-61)	2/cs