

PIPER SERIES.

The Piper series, the best selling walker in its class, was designed for those who need a lightweight yet solid and stable platform.

It's the one of the strongest aluminum walkers on the market today and is available in four sizes. Piper Super Mini with a 17" seat height weighing just 12 lbs to the Piper Tall with a 23" seat height weighing just 13 lbs.

FEATURES.

Aluminum Frame Superb stability, and style

Seat Solid plastic base with removable fitted waterproof cushion

Mesh Basket Removable

Wheels 8" 40 mm wide, Solid grip / Longer lasting polyurethane, Suitable for indoor & outdoor use

Colour Vibrant Red, Sapphire Blue, Silver, Black, Forest Green



SPECIFICATIONS.

Super Mini



12 lbs
16.5" Seat Width
17" Seat Height
21" Width x 26" Length
29" - 33" Handle Height
275 lbs Weight Capacity

Mini



13 lbs
18" Seat Width
19" Seat Height
21" Width x 26" Length
30" - 34" Handle Height
275 lbs Weight Capacity

Regular



13 lbs
18" Seat Width
21" Seat Height
21" Width x 26" Length
31" - 36" Handle Height
275 lbs Weight Capacity

Tall



13 lbs
18" Seat Width
23" Seat Height
23" Width x 26" Length
35" - 39" Handle Height
275 lbs Weight Capacity

OPTIONAL ACCESSORIES.

- Cup Holder
- Cane Holder
- Curb Climber
- Slow Down Brakes
- Deluxe Back Strap Pad
- Oxygen Tank Holder
- Mono Brakes
- Tote Bag (fits inside basket)
- Tray - Seat Mounted



Evolution
TECHNOLOGIES*

EVOLUTION TECHNOLOGIES INC.
2530 Davies Avenue,
Port Coquitlam, BC V3C2J9

tf +1 800.556.2558
info@evolution.ca
www.evolution.ca

EVOLUTION TECHNOLOGIES USA INC.
2657 Delta Ring Road
Ferndale WA 98248

tf +1 800.839.2558
info@evolutionwalker.com
www.evolutionwalker.com