

INVACARE®
TRACER® SX5
SERIES WHEELCHAIRS



The **Invacare® Tracer® SX5 Wheelchair's** lightweight frame, weighing less than 36 lb., makes this chair perfect for rental and long- or short-term use. This wheelchair offers multiple seat widths, arm styles, sectional back and seat-to-floor heights by utilizing the dual axle position frame base. The 20" and 22" widths include a semi-pneumatic caster to accommodate a weight capacity of 300 lb. Combining the design and technology of the Invacare® 9000 series wheelchairs make it possible to interchange components with the Tracer SX5 wheelchair.

The **Invacare® Tracer® SX5 Recliner Wheelchair** offers the durability of a lightweight manual wheelchair frame with the added versatility and comfort of a full reclining back. The Tracer SX5 Recliner wheelchair offers dynamic recline range from 90° to 180°. This chair is available with multiple seat widths, arm styles, front riggings, a full range of accessories and a 300 lb. weight capacity on 20" and 22" widths.

Invacare® Tracer® SX5 Wheelchair Features

Durable, low-maintenance, carbon steel frame is long-lasting

Dual axle positions allow variation of seat-to-floor heights

14-gauge crossbraces add strength and durability

Urethane rear tires mounted on "no flex" wheels offer superior performance

Sectional back available in 16" and 18" back heights

Low-maintenance, dual-sealed precision bearings

Heavy-duty inner liners keep seat and back from stretching

Meets RESNA standards for manual wheelchairs in product performance

Invacare® Tracer® SX5 Recliner Wheelchair Features

Dynamic range of recline from 90° to 180°

Durable, low-maintenance, carbon steel frame

Urethane rear tires, mounted on "no flex" wheels, offer superior performance

Low-maintenance, dual-sealed precision bearings

Heavy-duty inner liners keep seat and back from stretching

Upholstery is dual-embossed

Meets RESNA standards for manual wheelchairs in product performance



1. Dual axle positions
2. Urethane rear tires
3. Headrest
4. Seating options available
5. Builder options available



Invacare Corporation
www.invacare.com

USA

One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-8800

Canada

570 Matheson Blvd. E.,
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

Model #	Overall Width	Overall Height	Overall Depth w/o riggings	Folded Width	Seat Dimensions	Seat-to-Floor Height	Back Height
TRSX50, TRSX52 TRSX56, TRSX58	22.5"-30.5"	N/A	34"	N/A	14"-22"x16"	17.5", 19.5"	16", 18"
TRSX5RC6, TRSX5RC8	N/A	N/A	30.25"	N/A	14"-22"x16"	17.5", 19.5"	34"
	Caster Height	Wheel Height	Weight w/o Rigging	Weight Capacity	Shipping Weight	Carton Dimensions	Limited Warranty
TRSX50, TRSX52 TRSX56, TRSX58	6", 8"	20", 22", 24"	34 lb.	250 lb./300 lb.	43 lb.	32"x13"x 37"	5 Years
TRSX5RC6, TRSX5RC8	8"	20", 22", 24"	53 lb.	250 lb./300 lb.	62 lb.	36"x13"x 44"	5 Years

©2019

Invacare Corporation.
All rights reserved.
Trademarks are identified by the symbols ™ and ®.
All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Specifications are subject to change without notification.
Form No. 11-225
Rev. 04/14
131066

