




PARSONS

ADL



17S132

PARSONS MULTI PURPOSE POSITIONING BELT Ceinture de positionnement multifonction

This belt will fit a variety of bathlifts and hygiene chairs that offer bolt-on installation. Made of 2" (5cm) soft 100% nylon white webbing with a white plastic buckle. Attaches to any chair back utilizing the two aluminum adjustable belt anchor at each end of the belt. Undo a bolt on the back of the chair in order to attach the belt anchor. Maximum length approx 52" (132 cm).

Cette ceinture peut être attachée à une variété de leviers de bains et de chaises d'aisance sur lesquelles on peut la visser. Fabriquée de sangle douce blanche en 100 pour cent nylon 2 po (5 cm) de large, elle a une boucle fermeture rapide à pression en plastique blanche. Cette ceinture est vissée à l'arrière de la chaise en utilisant les deux supports d'ancrage en aluminium ajustables qui se trouvent sur les deux extrémités de la ceinture. Dévissez le boulon à l'arrière de la chaise pour y fixer le support d'ancrage de la ceinture. Longueur maximale est environ 52 po (132 cm).

Limited One Year Warranty

Parsons ADL warrants to the original owner of this product that it is free from defects in material and workmanship for one year from date of purchase. This warranty is offered provided the product has not been abused, misused, altered or damaged after purchase. (This includes damage due to the use of harsh chemicals.)

Garantie limitée d'un an

Parsons ADL garantit au propriétaire d'origine du présent produit que celui-ci est exempt de défauts de matériaux et de fabrication pour une durée d'une année à compter de la date d'achat. Cette garantie est offerte à condition que le produit n'ait pas été utilisé abusivement, mal utilisé, modifié ou endommagé après l'achat. (Ceci comprend les dommages résultant de l'utilisation de produits chimiques puissants.)

Parsons ADL Inc. 1986 Sideroad 15, Tottenham, Ontario L0G 1W0 Canada

Tel: 905.936.3580 - 1.800.263-1281 / Fax: 905.936.3585 - 1.888.225.5535

Email: custserv@parsonsadl.com

