

ACCESSORIES

CANE LIGHT



8-1041-00
TIP-SEE CLIP on LED light



Available in retail packaging

CANE STAYS



8-1039-00
Cane Clip
fits up to
1" diameter canes.

8-1040-08
Black Cane Stay,
fits up to 1" diameter canes.
8-1040-54
Lilac Cane Stay
fits up to 1" diameter canes.

HAND GRIPS



8-2199-08
Hypalon replacement grip, 3/4" inside diameter
by 12" long, black.

8-2197-08
Hypalon replacement grip, 7/8" inside diameter
by 13" long, black.



8-2198-08
Hypalon replacement grip,
3/4" inside diameter by
4" long, black.

8-2196-08
Hypalon replacement grip,
7/8" inside diameter by
5" long, black.



8-2195-08
Latex free, PVC replacement grip, 7/8" inside
diameter by 5" long, black.

8-2201-08
Plug & strap
for 3/4" tube.

8-2202-08
Plug & strap
for 7/8" tube.

ROPE WRIST STRAPS

Elastic loop. One size fits all.



	8-1001-08	Black
	8-1001-07	Brown
	8-1001-06	Blue
	8-1001-05	Hunter Green
	8-1001-04	Burgundy
	8-1001-03	Tan
	8-1001-02	Red
	8-1001-01	Purple

BRAIDED WRIST STRAPS

Elastic loop with clip. One size fits all.



	8-1015-08	Black
	8-1015-07	Brown

LEATHER STRAP



8-1025-08
Genuine
black
leather strap
for canes.



8-3203-80
Miter box and saw for sizing canes.
Comes with screwdriver tool.

9-0420-00
Measuring cane.
Aids in determining
patients cane
height without
numerous cutting
of non-adjustable
type canes. 28"-41"

**Note: not a weight
bearing cane**