





## **Shower Bench with Removable Padded Arms**

Item #	UOM
12440-2	2/cs
12440-1	1/cs

- Removable, padded armrests provide comfort and added security.
- Blow-molded bench provides comfort and strength.
- Drainage holes in bench reduce slipping.
- Aluminum frame is lightweight, durable and corrosion proof.
- · Leg height adjusts in 1" increments.

specifications	
Outside Legs	18.5" (W) x 18.25" (D)
Seat Dimensions	15.75" (W) x 13" (D)
Seat Height	16" - 21"
Weight (-1   -2)	5.5 lbs.   10.5 lbs.
Carton Shipping Weight (-1   -2)	8 lbs.   14 lbs.
Weight Limit	300 lbs.
Warranty	Limited Lifetime

## **Shower Stool**

Item #	<b>Product Description</b>	Color	UOM
RTL12004KD	Knocked Down	White	1/cs
12004KDRB-1	Knocked Down	Blue 📗	1/cs
12004KDRT-1	Knocked Down	Teal	1/cs

- Aluminum frame is lightweight, durable and corrosion proof.
- Impact-resistant, composite seat.
- Legs are height adjustable in 1" increments, and are crack proof and tarnish resistant.
- Easy, tool-free assembly.
- New pinch-free cover allows for push pins to be depressed without pinching fingers. The audible "click" will be heard when secure. (Figure A)
- Support collar prevents leg movement and reduces rattle. (Figure B)

specifications	
Outside Legs	12.5" (W) x 12.5" (D)
Seat Dimensions	13" (W) x 13" (D)
Seat Height	13.5" - 21"
Weight	3 lbs.
Carton Shipping Weight	5 lbs.
Weight Limit	300 lbs.
Warranty	Limited Lifetime